

## Mouse IL4 protein

Cat. No. GTX00401-pro

Species	Mouse
---------	-------

Package

20 µg

## 产品说明

## 摘要

Interleukin-4 Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 120 amino acids and having a molecular mass of 13.5 kDa.

## 属性

## 形式

Lyophilized powder

## 存储溶液

Batch dependent (please contact us for details)

## 存放说明

Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

## Region/Sequence

MHIHGCDKNH LREIIGILNE VTGEGTPCTE MDVPNVLTAT KNTTESELVC RASKVLRIFY LKHGKTPCLK KNSSVLMELQ RLFRAFRCCLD SSISCTMNES KSTSLKDFLE SLKSIMQMDY S.

## 表达系统

E. coli

## 纯度

&gt; 98% by SDS-PAGE.

## 内毒素

&lt; 1 EU/µg

## 偶联

Unconjugated

## 注意事项

仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).