

Mouse SCF protein

Cat. No. GTX00406-pro

Species	Mouse
---------	-------

Package	10 µg
---------	-------

产品说明

摘要

Stem Cell Factor Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 165 amino acids and having a molecular mass of 18.309 kDa.

属性

形式

Lyophilized powder

存储溶液

Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 10mM Sodium Citrate (pH4.0).

保存剂

No preservatives

存放说明

Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence

MKEICGNPVT DNVKDITKLV ANLPNDYMIT LNYVAGMDVL PSHCWL RDMV IQLSLSLTTL LDKFSNISEG LSNYSIIDKL GKIVDDL VLC MEENAPKNIK ESPKRPETRS FTPEEFSIF NRSIDAFKDF MVASDTSDCV LSSTLGPEKD SRVSVTKPFM LPPVA.

表达系统

E. coli

纯度

> 98% by SDS-PAGE and HPLC.

内毒素

< 1 EU/µg

偶联

Unconjugated

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).