

## CD8 alpha antibody [Hit8a] (APC)

Cat. No. GTX01486-07

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	FCM
种属反应	Human, Chimpanzee, Rhesus Monkey

引用文献 (3)

## 产品说明

## 摘要

The Hit8a antibody is specific for the 32-34 kDa alpha chain of human CD8, known as CD8a or CD8 alpha. CD8a can form a homodimer (CD8 alpha-alpha), but is more commonly expressed as a heterodimer with a second chain known as CD8b or CD8 beta. CD8 acts as a co-receptor for antigen recognition and subsequent T cell activation that is initiated upon binding of the T cell receptor (TCR) to antigen-bearing MHC Class I molecules. The cytoplasmic domains of CD8 provide binding sites for the tyrosine kinase lck, facilitating intracellular signaling events that lead to T cell activation, development, and cytotoxic effector functions. CD8+ cytotoxic T cells (CTLs) play an important role in inducing cell death of tumor cells, as well as cells infected by virus, bacteria or parasites. The Hit8a antibody is widely used as a phenotypic marker for CD8 on cytotoxic T cells, thymocytes, as well as on certain cell types that do not also express the TCR, including some NK cells and lymphoid dendritic cells. If used together with an alternative Anti-Human CD8a clone RPA-T8, the Hit8a antibody will not block binding of RPA-T8 to CD8a.

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

## Suggested dilution

## Recommended dilution

FCM	0.125 µg (5 µl) for 10 <sup>5</sup> -10 <sup>8</sup> cells in 100 µl sample per test
-----	--

以下为常规应用缩写的中文注解

WB: 免疫印迹  
ICC/IF: 细胞染色  
IHC-P: 石蜡切片  
IHC-Fr: 冰冻切片

## 属性

形式	Liquid
存储溶液	10mM NaH <sub>2</sub> PO <sub>4</sub> , 150mM NaCl, 0.1% Gelatin
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
浓度	0.025 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Allophycocyanin (APC) <a href="#">Wavelength</a>




For full product information, images and publications, please visit our [website](#).

Date 2026 / 04 / 07 Page 1 of 2

## 注意事项

仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

## 產品图片



## GTX01486-07 FCM Image

FACS analysis of human peripheral blood lymphocytes using GTX01486-07 CD8 alpha antibody [Hit8a] (APC).

Solid line : primary antibody

Dashed line : isotype control

antibody amount : 0.125  $\mu$ g (5  $\mu$ l)



For full product information, images and publications, please visit our [website](#).