

IL4 ELISA pair [IL4-I/IL4-II]

Cat. No. GTX03091

实验应用	ELISA, Sandwich ELISA
种属反应	Human

Package
1 pair

产品说明

Product Content

Cat No.	Type	宿主	偶联	Package
GTX02993	Capture antibody	Mouse	Unconjugated	250 µg
GTX02994-02	Detection antibody	Mouse	Biotin	250 µg

属性

形式 Liquid


存放说明

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

注意事项

仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX03091 ELISA Image

Sandwich ELISA analysis of human IL-4 protein using GTX02993 IL4 antibody [IL4-I] as coating antibody and GTX02994-02 IL4 antibody [IL4-II] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).