

ATP7b antibody [S62-29]

Cat. No. GTX03698

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ICC/IF, IP, IHC
种属反应	Human, Mouse, Rat

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:100
IP	Assay dependent
IHC	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


属性

形式	Liquid
存储溶液	PBS, 50% Glycerol
保存剂	0.09% sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片



GTXT03698 ICC/IF Image


ICC/IF analysis of 4%PFA-fixed NIH-3T3 using GTXT03698 ATP7b antibody [S62-29].

Green : Primary antibody

Blue : DAPI

Red : F-Actin

Dilution : 1:100



GTXT03698 ICC/IF Image


ICC/IF analysis of 4%PFA-fixed SH-SY5Y using GTXT03698 ATP7b antibody [S62-29].

Green : Primary antibody

Blue : Hoechst nuclear stain

Red : F-Actin

Dilution : 1:100



GTXT03698 WB Image

WB analysis of rat brain lysate using GTXT03698 ATP7b antibody [S62-29].

Dilution : 1:1000

Loading : 20 µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).