

Mouse Heart (Normal) tissue lysate

Cat. No. GTX03710

Species Mouse**Package**

100 µg

产品说明

摘要

Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

Tissue type

Normal

属性

形式

Liquid

存储溶液

HEPES, MgCl₂, KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, NP-40, and a cocktail of protease inhibitors

保存剂

No preservatives

存放说明

Store the protein lysate at -70°C for long term storage. After the protein is thawed, store it at 2 - 8°C for continuous use. Avoid freeze and thaw cycles.

浓度

5 mg/ml (Please refer to the vial label for the specific concentration.)

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).