

Mouse Anti-Human IgG4 (pFc') antibody [HP6023] (Biotin)

Cat. No. GTX04186-02

宿主	Mouse
克隆	Monoclonal
同种型	IgG3
实验应用	WB, IHC-Fr, FCM, IP, ELISA
种属反应	Human

Package
500 µg

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:5000-1:10000

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片


IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	PBS
保存剂	0.09% sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
免疫原	Human IgG4 myeloma protein
纯化方式	Purified IgG3
偶联	Biotin
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

產品圖片



GTX04186-02 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTX04186-02 Mouse Anti-Human IgG4 (pFc') antibody [HP6023] (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](#).