

## Goat Anti-Mouse IgG1 (Heavy chain) antibody, pre-adsorbed (FITC)

Cat. No. GTX04207-06

宿主	Goat
克隆	Polyclonal
同种型	IgG
实验应用	ELISA
种属反应	Mouse

Package  
1 mg

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ELISA	1:100-1:400

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


## 产品说明

Pre-adsorbed with Mouse IgG2a, IgG2b, IgG3, IgM, and IgA; human, bovine, and horse serum proteins (SP). May react with immunoglobulins from other species.

## 属性

形式	Liquid
存储溶液	PBS
保存剂	0.09% sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
免疫原	Pooled antisera from goats hyperimmunized with mouse IgG1
纯化方式	Purified by antigen-affinity chromatography
偶联	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

**產品图片**

**GTX04207-06 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04207-06 Goat Anti-Mouse IgG1 (Heavy chain) antibody, pre-adsorbed (FITC).



For full product information, images and publications, please visit our [website](#).