

Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (Biotin)

Cat. No. GTX04209-02

| | |
|------|------------|
| 宿主 | Goat |
| 克隆 | Polyclonal |
| 同种型 | IgG |
| 实验应用 | ELISA |
| 种属反应 | Mouse |

Package
1 mg

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ELISA | 1:5000-1:20000 |

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

产品说明


Reacts with the heavy chain of C57BL/6 mouse IgG2c. Pre-adsorbed with Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. May react with immunoglobulins from other species

属性

| | |
|------|-----------------------------------------------------------------------------------------|
| 形式 | Liquid |
| 存储溶液 | PBS |
| 保存剂 | 0.09% sodium azide |
| 存放说明 | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| 浓度 | 0.5 mg/ml (Please refer to the vial label for the specific concentration.) |
| 免疫原 | Pooled antisera from goats hyperimmunized with mouse IgG2c |
| 纯化方式 | Purified by antigen-affinity chromatography |
| 偶联 | Biotin |
| 注意事项 | 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。 |

For full product information, images and publications, please visit our [website](#).

產品圖片

**GTX04209-02 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04209-02 Goat Anti-Mouse IgG2c (Heavy chain) antibody, pre-adsorbed (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](#).