

## Goat Anti-Mouse IgG2a (Heavy chain) antibody, pre-adsorbed (AP)

Cat. No. GTX04213-03

宿主	Goat
克隆	Polyclonal
同种型	IgG
实验应用	ELISA
种属反应	Mouse

Package  
1 ml

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ELISA	1:2000-1:4000

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 产品说明


Pre-adsorbed with Mouse IgG1, IgG2b, IgG3, IgM, and IgA; human, bovine, and horse serum proteins (SP). May react with immunoglobulins from other species.

## 属性

形式	Liquid
存储溶液	50mM Tris, 1mM MgCl <sub>2</sub> , 50% Glycerol
保存剂	0.09% sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
免疫原	Pooled antisera from goats hyperimmunized with mouse IgG2a
纯化方式	Purified by antigen-affinity chromatography
偶联	Alkaline Phosphatase (AP)
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

產品圖片



**GTX04213-03 ELISA Image**

ELISA analysis of purified mouse immunoglobulins using serially diluted GTX04213-03 Goat Anti-Mouse IgG2a (Heavy chain) antibody, pre-adsorbed (AP).



For full product information, images and publications, please visit our [website](#).