

ATP5G1 antibody [N1C1]

Cat. No. GTX104671

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB
种属反应	Human

引用文献 (1)

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1:500-1:3000

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性


形式	Liquid
存储溶液	PBS, 40% Glycerol
保存剂	0.01% Thimerosal
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_1240492
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

Date 2025 / 09 / 09 Page 1 of 2

產品圖片



GTX104671 WB Image


Sample (30 µg whole cell lysate)

A:MOLT4 (GTX27912)

15% SDS PAGE

GTX104671 diluted at 1:1000

The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX104671 WB Image

The data was published in the journal Sci Rep in 2017. [PMID: 28701793](https://doi.org/10.1038/s41598-017-28701-7)



For full product information, images and publications, please visit our [website](https://www.genetex.com).