

PRPF31 antibody

Cat. No. GTX117081

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB, IHC-P
种属反应	Human, Mouse, Rat

引用文献 (1)

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


属性

形式	Liquid
存储溶液	PBS, 20% Glycerol
保存剂	0.01% Thimerosal
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_2885412
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

Date 2025 / 09 / 09 Page 1 of 2

產品圖片




GTX117081 IHC-P Image

PRPF31 antibody detects PRPF31 protein at nucleus speckles on rat fore brain by immunohistochemical analysis.

Sample: Paraffin-embedded rat fore brain.

PRPF31 antibody (GTX117081) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




GTX117081 IHC-P Image

PRPF31 antibody detects PRPF31 protein at nucleus speckles on mouse testis by immunohistochemical analysis.

Sample: Paraffin-embedded mouse testis.

PRPF31 antibody (GTX117081) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX117081 WB Image

PRPF31 antibody detects PRPF31 protein by western blot analysis.

A A549

B H1299

C HCT116 whole cell lysate/extract

7.5 % SDS-PAGE

PRPF31 antibody (GTX117081) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).