

ELK1 (phospho Ser383) antibody

Cat. No. GTX11848

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB
种属反应	Human, Mouse, Rat

实验应用

应用说明

WB: Use at a dilution of 1/1000. This concentration is determined using Elk1 fusion protein phosphorylated with purified Erk2. Low levels of endogenous Elk1 protein present in most cell types makes it difficult to detect by immunoblotting. Overexpression of Elk1 by transfection or immunoprecipitation, followed by immunoblotting will help improve detection. Predicted molecular weight: 45 kDa. Optimal dilutions/concentrations should be determined by the end user.

属性

形式	Liquid
存储溶液	PBS, 10mM HEPES, 150mM NaCl, 0.01% BSA, 50% Glycerol
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).