

PSD95 antibody [7E3-1B8]

Cat. No. GTX13552

宿主	Mouse
克隆	Monoclonal
同种型	IgG2a
实验应用	WB, ICC/IF, IP
种属反应	Mouse, Rat, Bovine

实验应用


应用说明

ICC: Use at an assay dependent dilution. IP: Use at an assay dependent dilution. WB: Use at a concentration of 1 µg/ml. Detects a band of approximately 100 kDa. Optimal dilutions/concentrations should be determined by the end user.

属性

形式	Liquid
存储溶液	PBS, 50% Glycerol
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

产品图片

**GTX13552 WB Image**

Western blot analysis of bovine brain tissue extract, probed with PSD-95, mAb (7E3-1B8).



For full product information, images and publications, please visit our [website](#).