

Rad9 antibody [93A535]

Cat. No. GTX13600

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ICC/IF, IHC-P, IP
种属反应	Human

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1 - 2 µg/ml
ICC/IF	1:10 - 1:500
IHC-P	1:10 - 1:500
IP	1:10 - 1:500

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片


IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	PBS, 0.05% BSA
保存剂	0.05% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).


產品圖片



GTX13600 WB Image

WB analysis of Rad9 recombinant protein using GTX13600 Rad9 antibody [93A535].

Dilution : 2 µg/ml



GTX13600 IHC-P Image

IHC-P analysis of human prostate adenocarcinoma tissue using GTX13600 Rad9 antibody [93A535].

Dilution : 5 µg/ml

Antigen retrieval : 10 mM sodium citrate buffer, pH 6.0



For full product information, images and publications, please visit our [website](https://www.genetex.com).