

Human Appendix (Normal) tissue lysate

Cat. No. GTX15355

Species Human**Package**

200 µg

产品说明

摘要

Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

Tissue type

Normal

属性

形式

Liquid

存储溶液

HEPES, MgCl₂, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

保存剂

No preservatives

存放说明

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度

5 mg/ml (Please refer to the vial label for the specific concentration.)

纯化方式

Whole Tissue Lysate

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).