

## Hepatitis C virus NS4A protein antibody [5D4/10E7]

## Cat. No. GTX19048

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, IHC-P, ELISA
种属反应	Hepatitis C virus

## 实验应用

## 应用说明

ELISA: Use at an assay dependent dilution. IHC-P: The dilution fold should be from 1:100 to 1:1000. WB: Use at a concentration of 0.5 µg/ml to detect 0.5 µg/lane of synthetic NS4 peptide and 0.1 µg/lane recombinant chimeric HCV polyprotein. Optimal dilutions/concentrations should be determined by the end user.

## 属性

形式	Liquid
存储溶液	PBS
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).