

Mouse Heart tissue lysate (14 days old)

Cat. No. GTX27194

Species Mouse**Package**

100 µg

产品说明

摘要

Heart tissue lysate (14 Days Old) was prepared by homogenization in modified RIPA buffer. Tissue and cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type

Normal

属性

形式

Liquid

存储溶液

1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

保存剂

No preservatives

存放说明

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度

2 mg/ml (Please refer to the vial label for the specific concentration.)

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).