

## TNF Receptor II antibody [HM102]

Cat. No. GTX27369

|      |  |
|------|--|
| 宿主   | Rat                                    |
| 克隆   | Monoclonal                             |
| 同种型  | IgG2a                                  |
| 实验应用 | WB, IHC-Fr, FCM, IP, ELISA, Activation |
| 种属反应 | Mouse                                  |

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

| Suggested dilution | Recommended dilution   |
|--------------------|--|
| WB                 | Assay dependent  |
| IHC-Fr             | Assay dependent  |
| FCM                | 10µg/ml antibody in PBS/1% serum / 10 <sup>5</sup> microglia cells |
| IP                 | Assay dependent  |
| ELISA              | Assay dependent  |
| Activation         | 2µg/ml   |

**Note : A reduced sample treatment and SDS-Page was used.  
In most cases acts agonistically.**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 属性

|      |   |
|------|---|
| 形式   | Liquid  |
| 存储溶液 | Filter-sterilized PBS, 0.1% BSA   |
| 保存剂  | No preservatives  |
| 存放说明 | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| 浓度   | 100 µg/ml (Please refer to the vial label for the specific concentration.)              |
| 偶联   | Unconjugated  |
| 注意事项 | 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。  |



For full product information, images and publications, please visit our [website](#).