

Rat Heart (Normal) tissue lysate

Cat. No. GTX27940

Species	Rat
---------	-----

Package	100 µg
---------	--------

产品说明

摘要

Heart tissue lysate was prepared by homogenization in modified RIPA buffer. Tissue and cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type

Normal

属性

形式

Liquid

存储溶液

1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

保存剂

No preservatives

存放说明

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度

2 mg/ml (Please refer to the vial label for the specific concentration.)

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).