

## BMF peptide

Cat. No. GTX29654

实验应用	Neutralizing /Inhibition
Species	Human

Package  
50 µg

## 实验应用

## 应用说明

The peptide is used for blocking the activity of anti-Bmf (C terminal) antibody (GTX29654). Incubating the peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the activity of anti-Bmf (C terminal) in Western blotting.

## 属性

形式	Liquid
存储溶液	PBS, 0.1% BSA
保存剂	0.02% Sodium azide
存放说明	Store the peptide at -20°C, stable for one year.
浓度	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).