

Mouse IL6 protein

Cat. No. GTX29731-pro

Species Mouse**Package**

10 µg

产品说明

摘要

Interleukin-6 Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 187 amino acids and having a molecular mass of 21.709 kDa.

属性

形式 Lyophilized powder**存储溶液** Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 0.2µm-filtered PBS (pH7.4).**保存剂** No preservatives**存放说明** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence**

FPTSQVRRGD FTEDTPNRP VYTTQVGGI ITHVLWEIVE MRKELCNGNS DCMNDDALA ENNLKLEIQ RNDGQYQYGY
NQEICLLKIS SGLLEYHSL EYMKNLNDN KDKARVLQR DTETLIHIFN QEVKDLHKIV LPTPISNALL TDKLESQKEW LRTKTIQFIL
KSLEEFKVT LRSTRQT.

表达系统 E. coli**纯度** > 98% by SDS-PAGE and HPLC.**内毒素** < 1 EU/µg**偶联** Unconjugated**注意事项**

For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).