

## Blocking Solution with Horse Serum (ready to use)

Cat. No. GTX30974

实验应用	WB, ICC/IF, ELISA, Blocking, IHC
Species	Horse

Package  
50 ml

## 产品说明

**摘要** Immuno Protein Blocking Solution with Horse Serum (ready to use)

## 实验应用

## 应用说明

The serum used to make this blocking solution is obtained from pools of healthy adult animals. The serum is heat inactivated at 56°C for 60 minutes and diluted and filtered. The serum is tested with appropriate biotinylated antibody for cross-reactivity. Immuno protein blocking solutions is ready to use and do not contain azide or phosphate buffers.

## 属性

**形式** Liquid**存放说明** Store at 2-8°C.**注意事项** 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).