

## CCDC141 blocking peptide

Cat. No. GTX31445-PEP

实验应用 Neutralizing /Inhibition

Package

50 µg

## 产品说明

**摘要** CCDC141 peptide is used for blocking the activity of CCDC141 antibody (catalog number GTX31445).

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

## Suggested dilution

Neutralizing /Inhibition

## Recommended dilution

Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 属性

**形式** Liquid**存储溶液** PBS, 0.1% BSA**保存剂** 0.02% Sodium azide**存放说明** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**浓度** 200 µg/ml (Please refer to the vial label for the specific concentration.)**偶联** Unconjugated**注意事项** 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。For full product information, images and publications, please visit our [website](#).