

AKAP 82 blocking peptide

Cat. No. GTX31595-PEP

实验应用 Neutralizing /Inhibition

Package
50 µg

产品说明

摘要 AKAP 82 peptide is used for blocking the activity of AKAP 82 antibody (catalog number GTX31595).

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution

Neutralizing /Inhibition

Recommended dilution

Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性

形式 Liquid

存储溶液 PBS, 0.1% BSA

保存剂 0.02% Sodium azide

存放说明 Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

浓度 200 µg/ml (Please refer to the vial label for the specific concentration.)

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).