

Human ATP6V1E1 protein

Cat. No. GTX48224-pro

Species	Human
---------	-------

Package

50 µg

属性

形式	Lyophilized powder
----	--------------------

存储溶液	Batch dependent (please contact us for details)
------	---

存放说明	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
------	---

Region/Sequence	a.a.1-226
-----------------	-----------

表达系统	E. coli
------	---------

纯度	.95
----	-----

偶联	Unconjugated
----	--------------

注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。
------	--



For full product information, images and publications, please visit our [website](#).