

## Mouse SCF protein (active)

Cat. No. GTX48296-pro

实验应用 Functional Assay

Species Mouse

Package

10 µg

## 实验应用

## 应用说明

The ED<sub>50</sub> was determined by the dose-dependent stimulation of the proliferation of human TF-1 cells is  $\leq 10.0$  ng/ml, corresponding to a specific activity of  $\geq 1 \times 10^7$  units/mg.

## 属性

形式 Lyophilized powder

存储溶液 Batch dependent (please contact us for details)

存放说明 Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MKEICGNPVT DNVKDITKLV ANLPNDYMIT LNYVAGMDVL PSHCWLRDMV IQLSLTLTLDKFSNISEG LSNYSIIDKL GKIVDDLVL MEENAPKNIK ESPKRPETRS FTPEEFSIF NRSIDAFKDF MVASDTSDCV LSSTLGPEKD SRVSVTKPFM LPPVA

表达系统 E. coli

纯度 &gt; 98% by SDS-PAGE and HPLC.

内毒素 &lt; 1 EU/µg

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).