

Rabbit Anti-Human Ig light chain antibody [RM129] (Biotin)

Cat. No. GTX60920

宿主	Rabbit
克隆	Monoclonal
同种型	IgG
实验应用	ICC/IF, IHC-P, FCM, ELISA, Sandwich ELISA
种属反应	Human

Package

50 µg

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ICC/IF	0.5µg/mL - 2µg/mL
IHC-P	0.5µg/mL - 2µg/mL
FCM	Assay dependent
ELISA	0.02µg/mL - 0.25µg/mL
Sandwich ELISA	Assay dependent

**Note : Capture : GTX20291/GTX60861/GTX60862/GTX60863/GTX60864/GTX60865
/GTX60866/GTX60867/GTX60868/GTX60869/GTX60870/GTX60873,
Detection : GTX60920**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

产品说明

This antibody reacts to both kappa and lambda light chains of human immunoglobulins. The antibody does not react to monkey (Cyno or Rhesus) IgG, mouse IgG, rat IgG, or goat IgG.

属性





形式	Liquid
存储溶液	PBS, 1% BSA, 50% Glycerol
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)



For full product information, images and publications, please visit our [website](#).

Date 2025 / 09 / 06 Page 1 of 2

免疫原	Human IgG
纯化方式	Protein A purified From tissue culture supernatant
偶联	Biotin
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

產品图片


For full product information, images and publications, please visit our [website](https://www.genetex.com).