

Respiratory Syncytial virus type B G protein antibody [HL1733]

Cat. No. GTX637381

宿主	Rabbit
克隆	Monoclonal
同种型	IgG
实验应用	WB, ICC/IF, IHC-P (cell pellet)
种属反应	Respiratory syncytial virus type B

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	Assay dependent
IHC-P (cell pellet)	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


属性

形式	Liquid
存储溶液	PBS
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片


**GTX637381 ICC/IF Image**

Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] detects Respiratory Syncytial virus type B Glycoprotein protein by immunofluorescent analysis.


Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: Respiratory Syncytial virus type B Glycoprotein stained by Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] (GTX637381) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTX30920).

**GTX637381 WB Image**

Non-transfected (-) and transfected (+) 293T whole cell extracts (30 μg) were separated by 12% SDS-PAGE, and the membrane was blotted with Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] (GTX637381) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.

**GTX637381 IHC-P (cell pellet) Image**

Respiratory Syncytial virus type B G protein antibody [HL1733] detects Respiratory Syncytial virus type B G protein protein by immunohistochemical analysis.

Sample: Paraffin-embedded mock and RSV Type B Major surface glycoprotein G transfected 293T.

Green: Respiratory Syncytial virus type B G protein stained by Respiratory Syncytial virus type B G protein antibody [HL1733] (GTX637381) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTX30920).

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](#).