

Taf4 (S. cerevisiae) antibody

Cat. No. GTX64179

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB
种属反应	Saccharomyces cerevisiae

实验应用

应用说明

Western blotting. (1:1,000-1:5,000 dilution)

属性

形式	Liquid
存储溶液	Serum
保存剂	0.1% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX64179 WB Image

Detection of Taf4p by Western blotting Lane 1: Extract of budding yeast The antiserum was diluted 5,000 fold before use. The molecular mass of Taf4 is 48 kDa



For full product information, images and publications, please visit our [website](#).