

## Varicella Zoster virus IE62 antibody [62B]

Cat. No. GTX64190

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ICC/IF, IP
种属反应	Varicella Zoster virus

引用文献 (1)

## 实验应用


## 应用说明

Western blotting (1:2,000-1:5,000), Immunoprecipitation (1:100), Immunofluorescence staining and Immunocytochemistry (1:50-1:100)

## 属性

形式	Liquid
存储溶液	Filter-sterilized PBS, 50% Glycerol
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

## 产品图片




## GTX64190 ICC/IF Image

Immunofluorescence staining of VZV IE62 protein in VZV-infected MRC-5 cells by using anti-VZV IE62 antibody (clone 62B). Anti-VZV IE62 antibody was used at 1/100 dilution. As second antibody, Alexa Fluor 488 donkey anti-mouse IgG [H+L] was used at 1/200 dilution. Nuclei were stained with Hoechst 33342.



For full product information, images and publications, please visit our [website](#).



#### GTX64190 WB Image

Identification of IE62 protein in VZV-infected cells by western blotting using anti-VZV IE62 antibody (clone 62B). Lane 1; VZV strain pOka infected MRC-5 cell lysate Lane 2; MRC-5 cell lysate ( uninfected negative control ) The anti-VZV IE62 antibody was used at 1/5,000 dilution



For full product information, images and publications, please visit our [website](#).