

Mouse TIM-4 protein, mouse IgG2a Fc tag (active)

Cat. No. GTX65777-pro**实验应用** Functional Assay**Species** Mouse**Package**

100 µg

实验应用**应用说明**

Measured by its ability to inhibit anti-CD3-induced proliferation of stimulated human T cells.

属性**形式** Lyophilized powder**存储溶液** Reconstitute with 100 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.**保存剂** No preservatives**存放说明** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of mouse Tim-4 (a.a. 22-279) is fused to the N-terminus of the Fc region of mouse IgG2a.**表达系统** HEK293 cells**纯度** ??98% by SDS-PAGE.**内毒素** < 5 EU/mg**偶联** Unconjugated**注意事项** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).