

Human NMT1 protein, His tag

Cat. No. GTX67597-pro

Species Human

Package
100 µg

实验应用

Observed MW 57 kDa.

属性

形式 Liquid

存储溶液 PBS

保存剂 No preservatives

存放说明 Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSSLVPRGSH MGSMADSEET AVKPPAPPLP QMMEGNGNGH
 EHCSDCENEE DNSYNRGGLS PANDTGAKKK KKKQKKKKEK GSETDSAQDQ PVKMNSLPAE RIQEIQAIE LFSVGQGPAK
 TMEEASKRSY QFWDTQPVPK LGEVNTGHP VEPDKDNIRQ EPYTLPQGFT WDALDLGDRG VLKELYTLN ENYVEDDDNM
 FRFDYSPEFL LWALRPPGWL PQWHGCVRVV SSRKLVGFIS AIPANIHIDYD TEKKMVEINF LCVHKKLRSK RVAPVLIREI TRRVHLEGIF
 QAVYTAGWL PKPVGTCRYW HRSLNPRKLI EVKFSHLRN MTMQRTMKLY RLPETPKTAG LRPMETKDIP VHQLLTRYL
 KQFHLPVMS QEEVEHWFYP QENIIDTFV ENANGVTFDF LSFYTLPTI MNHPTHKSLK AAYSFYNVHT QTPLLDLMSD
 ALVLAKMKGK DVFNALDLME NKTFLKLF GIGDGNLQYY LYNWKCPMSG AEKVGLVLQ

表达系统 E. coli

纯度 > 95% by SDS-PAGE.


偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片

**GTX67597-pro Image**

3 μ g Human NMT1 protein (GTX67597-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).