

## RAD54L antibody [38.4]

**Cat. No. GTX70305**

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ELISA
种属反应	Human

## 实验应用


## 应用说明

Optimal dilutions/concentrations should be determined by the researcher.

## 属性

形式	Liquid
存储溶液	PBS
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_372860
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

## 产品图片



## GTX70305 WB Image

Detection of human Rad54 protein using GeneTex Rad54 38.4 monoclonal antibody (GTX70305) in HeLa whole cell extract.



For full product information, images and publications, please visit our [website](#).