

Rat IgG2a isotype control (Biotin)

Cat. No. GTX74766

宿主	Rat
克隆	Monoclonal
同种型	IgG2a
实验应用	FCM

Package
100 µg

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
FCM	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

产品说明

This antibody is suitable for use as a negative control to assess non-specific binding of Rat IgG2a antibodies to mouse and dog cells. And It is not suitable for use as a negative control on human.

属性

形式	Liquid
存储溶液	PBS, 1% BSA
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
纯化方式	Protein A purified
偶联	Biotin
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).