

Fas antibody [LOB 3/17]

Cat. No. GTX75106

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	FCM, IP
种属反应	Human

实验应用


应用说明

For flow cytometry: Use 10ul of the suggested working dilution (1:50-1:100) to label 10^6 cells in 100ul. For IP: Use at 20 μ g/ml. Optimal dilutions/concentrations should be determined by the researcher.

属性

形式	Liquid
存储溶液	Tris, NaCl
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_376599
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX75106 FCM Image

Staining of HUT78 T cell line with Mouse anti Human CD95:FITC (GTX75107)



For full product information, images and publications, please visit our [website](#).