

Pseudomonas aeruginosa serotype 5c antibody [SD6930]

Cat. No. GTX80347

| | |
|------|------------------------|
| 宿主 | Mouse |
| 克隆 | Monoclonal |
| 同种型 | IgG1 |
| 实验应用 | ICC/IF, EIA |
| 种属反应 | Psuedomonas aeruginosa |

实验应用

应用说明

Optimal dilutions should be determined experimentally by the researcher.

属性

| | |
|------|--|
| 形式 | Liquid |
| 存储溶液 | PBS |
| 保存剂 | 0.1% Sodium azide |
| 存放说明 | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| 浓度 | Batch dependent (Please refer to the vial label for the specific concentration.) |
| 偶联 | Unconjugated |
| 注意事项 | 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。 |



For full product information, images and publications, please visit our [website](#).